

Hot plate DPH 40


Hot plate for slides



Hot plate for drying tissue sections on slides

COD. 53DPH402025

Manufacturer: Diapath

**Description**

DPH 40 hot plate is used into histology, pathology, cytology and bacteriology labs, for drying tissue sections slides.

Technical features

Overall dimensions:	250x200x80 mm (WxDxH)
Hot plate dimensions:	250x200 mm (WxD)
Overall weight:	3 Kg
Certification:	IVD – CE
Power:	230V – 50/60 Hz single-phase
Electrical power consumption:	100 Watt
Working temperature:	from 30°C to 70°C
Operational temperature:	from 50°C to 40°C
Humidity operational level:	60-65% max